

BRL R 1312

# BRL

AD635655

REPORT NO. 1312

## APPLICATION OF FAST NEUTRON ACTIVATION ANALYSIS TO POLYMER DOSIMETRY

by

Arthur D. Coates  
Joe A. Swisher

March 1966

Distribution of this document is unlimited.

U. S. ARMY MATERIEL COMMAND  
**BALLISTIC RESEARCH LABORATORIES**  
ABERDEEN PROVING GROUND, MARYLAND

Destroy this report when it is no longer needed.  
Do not return it to the originator.

The findings in this report are not to be construed as an official Department of the Army position, unless so designated by other authorized documents.

The use of trade names or manufacturers' names in this report does not constitute indorsement of any commercial product.

BALLISTIC RESEARCH LABORATORIES

REPORT NO. 1312

MARCH 1966

Distribution of this document is unlimited.

APPLICATION OF FAST NEUTRON ACTIVATION ANALYSIS TO POLYMER DOSIMETRY

Arthur D. Coates  
Joe A. Swisher

Terminal Ballistics Laboratory

PROPERTY OF U.S. ARMY  
GPO: 1966-000000  
BRL, INC., D. 21006

RDTE Project No. 1N022601A086

ABERDEEN PROVING GROUND, MARYLAND

B A L L I S T I C   R E S E A R C H   L A B O R A T O R I E S

REPORT NO. 1312

ADCoates/JASwisher/jrl  
Aberdeen Proving Ground, Md.  
March 1966

APPLICATION OF FAST NEUTRON ACTIVATION ANALYSIS TO POLYMER DOSIMETRY

ABSTRACT

Fast neutron activation analysis was used to measure the gamma dose absorbed by a polymer. A series of low density polyethylene films was irradiated by cobalt-60 gamma rays in the presence of air. The gamma doses ranged from  $1 \times 10^5$  to  $4.7 \times 10^7$  rads. The oxygen content of the polymer films increased with increasing dose. The increase in oxygen content of each film was monitored at  $1720 \text{ cm}^{-1}$  by infrared spectroscopy. The oxygen increase also was monitored by measuring the activation of oxygen-16. This latter method employs fast neutron activation of oxygen-16 and the measurement of the emitted radiation to determine gamma dose. The results of both methods of dosimetry are presented, and the improvements offered by the neutron activation technique are indicated.

Page intentionally blank

Page intentionally blank

Page intentionally blank

## TABLE OF CONTENTS

	Page
ABSTRACT. . . . .	3
INTRODUCTION. . . . .	7
EXPERIMENTAL PROCEDURES AND APPARATUS . . . . .	11
RESULTS AND DISCUSSION. . . . .	13
CONCLUSIONS . . . . .	17
REFERENCES. . . . .	19
DISTRIBUTION LIST . . . . .	21

Page intentionally blank

Page intentionally blank

Page intentionally blank

## INTRODUCTION

This report presents the results of a study on the application of fast neutron activation analysis to high level gamma ray dosimetry. The need for high level dosimetry has become greater with the advent of increasingly larger irradiation devices which emit gamma doses greater than  $1 \times 10^6$  rads in less than one hour. Previous studies by the Ballistic Research Laboratories utilized the radiation induced chemical changes in polyethylene to determine gamma ray dose by infrared spectral techniques.<sup>1\*</sup> This report discusses a method for extending the range of polymer dosimetry above  $5 \times 10^7$  rads and for improving the accuracy at lower dose levels. The method employs fast neutron activation of oxygen-16 and the measurement of radiation emitted during the decay processes for determining the gamma ray dose. Both methods of dosimetry are based upon the finding that the oxygen content of the low density polyethylene, when irradiated in air, increases with gamma ray dose.

Several chemical reactions occur within polyethylene when this material is exposed in air to gamma radiation.<sup>2</sup> These reactions include crosslinking, formation of random double bonds, hydrogen evolution, main and side chain cleavage and degradative oxidation.

Oxygen reacts with polyethylene when the latter is undergoing gamma irradiation, and forms several carbonyl and hydroxyl groups. The increase in oxygen content of polyethylene is observed spectrally at three different wavelength regions. The spectral regions are listed below:

3300 - 3400 $\text{cm}^{-1}$	Hydroxyl Groups
1600 - 1800 $\text{cm}^{-1}$	Carbonyl Groups
800 - 1200 $\text{cm}^{-1}$	Fingerprint Region

As the oxygen content of polyethylene increases, the spectral absorption increases. This change in absorption has been correlated with dose and is the basis of the technique for utilizing polyethylene film as a gamma ray detector. The use of infrared spectral techniques is a convenient and reasonably sensitive method for observing the increase in oxygen content

---


\* Superscript numbers denote references which may be found on page 19.




as a function of gamma dose. Figure 1 shows the change in spectral absorption at  $1720\text{ cm}^{-1}$  in gamma irradiated polyethylene versus gamma dose.

Published reports<sup>3-6</sup> have shown that fast neutron activation is a rapid and sensitive method for determining oxygen-16. The minimum detectable limits are reported to be 10 parts per million (ppm) in samples containing 10 mg or more of oxygen.<sup>4</sup> The extension of fast neutron activation analysis to polymer dosimetry, by utilizing the oxygen-16 contained in irradiated hydrocarbon polymers, should reduce the time required for determination of dose, particularly where large numbers of individual measurements are to be performed. In addition, it may be possible to extend the upper and lower limits of the useful range, and conceivably to improve the accuracy of the polymer dosimetry technique. The threshold energy required for neutron induced activation of oxygen-16 is 9.6 MeV.<sup>3</sup> The abundance of oxygen-16 in oxygen is 99.76 percent.<sup>7</sup>

It is convenient and reasonably inexpensive to obtain 14.3 MeV neutrons by using a Cockcroft-Walton type accelerator to accelerate deuterons to approximately 150 KeV. The deuterons then strike a tritium target.<sup>8</sup> The activation studies reported herein were performed with a system of this type. The fast neutrons are produced by the reaction:



The concentration of oxygen-16 in the irradiated polymer was determined from the quantity of radiation emitted by the decay of nitrogen-16. The half life of nitrogen-16 is 7.4 seconds. The reactions that occur during the activation of oxygen-16 are:



The disintegration scheme for nitrogen-16 is shown in Figure 2.




FIGURE 1. THE CHANGE IN ABSORBANCE AT  $1720 \text{ CM}^{-1}$  IN POLYETHYLENE VS. GAMMA DOSE.




FIGURE 2 — NITROGEN-16 DISINTEGRATION SCHEME.<sup>9</sup>

Reaction 3 is the basis for the analytical determination, since the quantity of the emitted gamma radiation is proportional to the oxygen concentration. This gamma radiation is detected and measured by scintillation spectrometry. The normalized activity data obtained from each activated polymer sample is then related to the known gamma dose absorbed previously by the respective polymer films in a high intensity gamma irradiation unit. The normalized counts are then plotted versus gamma dose; these provide the instrument calibration in counts as a function of gamma dose.

The specific oxygen content of each irradiated polymer film may be determined if desired. This is accomplished by relating the normalized counts from each activated sample with the data obtained in an identical manner from standard samples of known oxygen content using the relationship:

$$WT_x = \frac{A_x \times WT_o}{A_o} \quad (4)$$

where:

$A_o$  = total counts from oxygen standard

$A_x$  = total counts from unknown sample

$WT_o$  = weight of oxygen in standard sample

$WT_x$  = weight of oxygen in unknown sample

The weight of oxygen in each irradiated sample may be determined, and related to gamma dose received by the polymer film. This provides a method for determining gamma ray dose and has potential use in inter-laboratory comparisons where different types of radiation counting equipment are employed. The use of the oxygen weight characteristic also would permit laboratories without activation facilities to conduct polymer dosimetry studies by employing techniques such as the classical Unterzaucher method for oxygen determination.

It should be noted, that in the use of low density polyethylene for gamma dosimetry, each lot of polyethylene studied by BRL has required gamma dose calibration. Therefore, the oxygen content data should be considered typical only for a specific type and lot of polyethylene film.

#### EXPERIMENTAL PROCEDURES AND APPARATUS

The polyethylene film used in this dosimetry study was supplied by the E. I. DuPont De Nemours and Company on special request. The film had a uniform thickness of  $2.0 \pm 0.2$  mils and a density range between 0.915

and  $0.925 \text{ gm/cm}^3$ . The samples were prepared in strips 1.9 cm wide by 122 cm in length and weighed prior to irradiation. The infrared spectrum of each polymer film also was obtained prior to irradiation. The sample films were then placed in an open container, irradiated and subsequently reweighed.

The gamma irradiations were performed in the BRL high intensity gamma irradiator, which is a version of the AECL Gamma-cell system. The gamma dose in the irradiation cavity of this cobalt-60 unit was determined by the manufacturer and verified by these Laboratories. The calibrations and verifications were performed using the Fricke chemical dosimeter.<sup>10</sup> Isodose curves for each irradiation position in the cavity were established. The dose rate at the center of the cavity during the polymer exposures was  $1.08 \times 10^6 \text{ rads/hr} \pm 3.5 \text{ percent}$ .

Following the irradiation, an infrared spectrum was obtained of each polyethylene film. The change in spectral absorption at  $1720 \text{ cm}^{-1}$  also was obtained in order to monitor the oxygen increase attributable to carbonyl linkages. These measurements were obtained using the Perkin-Elmer Corporation Model 13 double beam, ratio recording spectrophotometer equipped with sodium chloride optics and densitometer attachment.

Each polymer film was then wrapped as a tight concentric unit and placed in an activation sample container in preparation for the fast neutron activation analysis. The polyethylene sample containers have been described elsewhere.<sup>8</sup>

The activation analysis of each of the gamma irradiated films was performed at the Oak Ridge National Laboratories. A Texas Nuclear Corporation Model 9501, 150 KV Cockcroft-Walton type generator was used. The beam current was 500  $\mu\text{a}$  and the target was tritiated titanium. A Radiation Instrument Development Laboratories 400 channel analyzer with an electric typewriter printout system was used in conjunction with two  $3" \times 3"$  NaI(Tl) scintillation crystals, a rapid transfer system for the activated sample and a lead shielded sample counting chamber.<sup>8</sup>

In these experiments the fast neutron bombardment of the sample had a duration of ten seconds, the transfer time with preset counting delay was three seconds and the counts from the sample were accumulated for thirty seconds which is approximately four half-lives of nitrogen-16. The lower and upper level discriminators of the multichannel analyzer were adjusted for optimum counting of the high energy gammas, and minimizing the counts from interfering reactions.

In order to compare gamma dose with the specific oxygen-16 content of the polymer films, samples of silicon dioxide were used as oxygen standards. The oxygen content of each irradiated film was determined by comparison of the activity of the gamma irradiated film with the activity of the oxygen standards.

## RESULTS AND DISCUSSION

The activity data obtained using fast neutron activation analysis in polymer dosimetry are presented in Figure 3. The neutron induced activity of each polymer film is plotted as normalized counts obtained from each sample versus gamma dose. The activity measured in each sample is attributable primarily to the radioactive decay of nitrogen-16 to stable oxygen-16. The quantity of radiation detected represents the total integrated counts obtained over the range of 1.7 to 9.7 MeV, which includes the 6.13 and 7.13 MeV gamma photopeaks of nitrogen-16. Using the techniques described herein, no significant interfering reactions were found in the neutron activated polymer films. Each activated film decayed following the half-life pattern of  $^{16}_7\text{N}$ .

The results show that the greater the gamma dose received by the polymer film, the larger the total number of counts detected subsequently in the neutron activated polymer film. This indicates that the oxygen-16 content of the polymer films increases with gamma dose, since it has been established that the number of counts obtained from the activated oxygen compounds are proportional to oxygen concentration. It is known that polyethylene irradiated in air by gamma rays contains radiation induced quantities of oxygen.<sup>11</sup> Figure 4 shows the percentage of oxygen present in irradiated polyethylene as determined by neutron activation.




FIGURE 3. THE COUNTS OBTAINED FROM ACTIVATED POLYETHYLENE SAMPLES VS. GAMMA DOSE.




FIGURE 4. RADIATION INDUCED OXYGEN CONTENT OF POLYETHYLENE VS. GAMMA DOSE.



Each gamma irradiated polymer film was studied by infrared analysis prior to neutron activation. Figure 1 presented earlier shows the increase in absorption at  $1720\text{ cm}^{-1}$  versus gamma dose, this change in absorption is attributable to radiation induced C=O formation in the polyethylene. The technique for relating dose with infrared absorption was reported previously, and is included here for comparison.

Polymeric films of the olefin type offer several advantages for use in dosimetry studies in the megarad range. Low density polyethylene in particular, has been found to be useful because of its relatively high degree of stability, minimum significance of post irradiation reactions, and simplified infrared spectrum. Dosimetry measurements with polyethylene films utilizing either activation analysis or infrared spectral techniques possess several desirable characteristics. The advantages include: (1) the time required for dose determinations by either method is less than five minutes; (2) the irradiated polymer detector may be analyzed at a given time and verified if desired at a later date; (3) the gamma dose to the detector may be provided in either single or multiple exposures; (4) the use of lengthy analytical procedures and highly purified compounds is not necessary and (5) the data obtained are relatively simple to resolve by utilizing predetermined calibration curves. In addition the data indicate that the polymer neutron activation technique has the added advantage of increased accuracy and extended dose range.

Figures 1 and 4 show respectively the infrared and activation dosimetry data using the same irradiated polymer film for both determinations. Over the dose range studied, it is observed that the experimental points on the activation plot show less random deviation than those of the infrared curve.

With either technique, the polymer films of the type used in these experiments have a lower usable dose limit in the range of  $1 \times 10^5$  rads. This limit results from sensitivity characteristics of the instrumentation, and not from the inherent characteristics of the polymer films. The data shown in Figure 3 indicate that the range may be extended by using activation analysis.

With respect to higher dose level determinations, the infrared technique is limited to measurable optical density changes. At  $1720\text{ cm}^{-1}$  the polyethylene film used here approaches an optical density of infinity at dose levels of approximately  $1 \times 10^8$  rads. Polyethylene film exposed to this dose level has begun to deteriorate physically and has become brittle. Thus, the determination of absorbed dose becomes increasingly difficult. Therefore, the analytical sensitivity and the physical condition of the polymer limit the upper usable dose range of the infrared technique.

Fast neutron analysis of the higher dose gamma irradiated polyethylene film is not limited significantly by instrument sensitivity. The accuracy of dose determinations has been found to increase with dose, and is approximately  $\pm 5$  percent at  $4.7 \times 10^7$  rads. It is believed that the upper gamma dose level for activation dosimetry will extend above  $5 \times 10^8$  rads. This assumption is based upon the recent finding that the oxygen concentration in irradiated film has been found to increase with dose up to  $5 \times 10^8$  rads, which was the largest dose studied. These latter analyses were performed using the classical Unterzaucher technique. Deterioration of polyethylene film caused by gamma radiation does occur but does not restrict the use of polymer films for gamma dosimetry, when using the activation technique.

### CONCLUSIONS

This report describes a gamma dose measurement technique which uses a polymer detector. The technique demonstrates the feasibility of utilizing fast neutron activation analysis as a tool for determining the gamma dose received by a polyethylene film. This technique requires that the polymer detector be irradiated in air. The oxygen uptake of the polymer detectors, by chemical reaction, is induced by ionizing radiation and is related to gamma dose. The determination of dose is performed subsequently, by correlation of the fast neutron induced activity of oxygen in the polyethylene films, with previously established calibrated gamma doses. The results have shown that the oxygen-16 activation technique used in conjunction with the polymer detector offers several advantages for gamma dosimetry use in the megarad region.

In addition it is shown that the oxygen content of low density polyethylene film can be determined by fast neutron activation analysis. It should be noted that the technique should in no way be limited to use with gamma irradiated polyethylene, but may be used to determine the oxygen content of any polyethylene regardless of the origin of the contained oxygen. It is also evident that the polyethylene films did react with sufficient atmospheric oxygen to be analyzed quantitatively for oxygen by neutron activation. The infrared spectra of the irradiated films show that the radiation induced oxygen content is present primarily in the form of carbonyl and hydroxyl groups.

As the course of the activation experiments progressed, the results indicated the potential usefulness of the technique for high level dosimetry. It became evident also, that other areas should be studied and some additional aspects amplified. Thus, it would be desirable to: (1) establish the optimum geometry and mass of the polymer detectors for maximum gamma sensitivity and widest useful operating range; (2) establish the overall accuracy of the optimized system; (3) investigate the gamma energy and rate dependence of the system; (4) establish the sensitivity of the system to other types of ionizing radiation; (5) investigate the response of the polymer detectors in pulsed and steady state mixed radiation fields, and (6) study polypropylene as an air irradiated gamma detector film for use with activation analysis.

#### ACKNOWLEDGEMENT

The authors wish to express their appreciation to Dr. J. E. Strain, Oak Ridge National Laboratories and Dr. H. E. Banta, Oak Ridge Institute of Nuclear Studies for their valuable technical assistance and guidance. The authors wish to acknowledge the experimental assistance given by W. S. Hayden of the Ballistic Research Laboratories.

ARTHUR D. COATES

JOE A. SWISHER

#### REFERENCES

1. Swisher, J. A. and Coates, A. D., A Study of the Dosimetric Properties of Low Density Polyethylene, Ballistic Research Laboratories Memorandum Report No. 1336, March 1961.
2. Black, R. M. and Charlesby, A., International Journal of Applied Radiation and Isotopes, 7, 126 (1959).
3. Veal, D. J. and Cook, C. F., Analytical Chemistry, 34, 178 (1962).
4. Steele, E. L. and Meinke, W. W., Analytical Chemistry, 34, 185 (1962).
5. Volborth, A. and Banta, H. E., Analytical Chemistry, 35, 2203 (1963).
6. Anders, O. U. and Birden, D. W., Analytical Chemistry, 36, 287 (1964).
7. Handbook of Chemistry and Physics, Forty-Fourth Edition, Chem. Rubber Publishing Co., Cleveland, Ohio (1962).
8. Strain, J. E., Hampton, W. J. and Liddicotte, G. W., The ORNL Analytical Chemistry Division 150-KV Cockcroft-Walton Generator, ORNL-TM 362, September 1962.
9. Radiological Health Handbook, PB 12174R, U. S. Dept. of Health, Education and Welfare, Revised Edition, September 1960.
10. International Conference on Peaceful Uses of Atomic Energy, 14, 179 (1955).
11. Charlesby, A., Atomic Radiation and Polymers, Pergamon Press, New York, N.Y. (1960).

Page intentionally blank

Page intentionally blank

Page intentionally blank

# DISTRIBUTION LIST

<u>No. of</u> <u>Copies</u>	<u>Organization</u>	<u>No. of</u> <u>Copies</u>	<u>Organization</u>
20	Commander Defense Documentation Center ATTN: TIPCR Cameron Station Alexandria, Virginia 22314	1	Commanding General U.S. Army Mobility Command Warren, Michigan 48090
1	Director of Defense Research & Engineering (OSD) ATTN: Asst Dir/Tac Msl & Ord Washington, D.C. 20301	1	Commanding General U.S. Army Tank-Automotive Center ATTN: Mr. J. G. Brookes Warren, Michigan 48090
3	Director Defense Atomic Support Agency ATTN: Col C. Daniels Washington, D.C. 20301	1	Commanding Officer U.S. Army Aviation Materiel Laboratories Fort Eustis, Virginia 23604
1	Commander, Field Command Defense Atomic Support Agency ATTN: FC Tech Lib Sandia Base P.O. Box 5100 Albuquerque, New Mexico 87115	1	Commanding Officer U.S. Army Engineer Research & Development Laboratories ATTN: STINFO Div Fort Belvoir, Virginia 22060
1	Commanding General U.S. Army Materiel Command ATTN: Dir of Rsch & Labs Washington, D.C. 20315	1	Commanding General U.S. Army Munitions Command Dover, New Jersey 07801
1	Commanding General U.S. Army Materiel Command ATTN: AMCRD-RP-B Washington, D.C. 20315	1	Commanding Officer U.S. Army Frankford Arsenal ATTN: Pitman-Dunn Lab Lib Br, 0270, Bldg 40 Philadelphia, Pennsylvania 19137
1	Commanding General U.S. Army Electronics Command ATTN: AMSEL-RD-H Fort Monmouth, New Jersey 07703	1	Commanding Officer U.S. Army Picatinny Arsenal Dover, New Jersey 07801
1	Commanding General U.S. Army Missile Command ATTN: Mr. C. G. Morgan Redstone Arsenal, Alabama 35809	1	Commanding Officer U.S. Army Dugway Proving Ground ATTN: Radl Op Div Dugway, Utah 84022
		1	Commanding General U.S. Army Weapons Command ATTN: SWERI-RDL, Mr. Eisler Rock Island, Illinois 61202

# DISTRIBUTION LIST

<u>No. of</u> <u>Copies</u>	<u>Organization</u>	<u>No. of</u> <u>Copies</u>	<u>Organization</u>
1	Commanding Officer U.S. Army Materials Research Agency ATTN: Dr. L. S. Foster Watertown, Massachusetts 02172	1	Commanding Officer U.S. Army Research Office ATTN: Mr. R. Mace Box CM, Duke Station Durham, North Carolina 27706
1	Commanding Officer U.S. Army Nuclear Defense Laboratory ATTN: Dr. H. Donnert Edgewood Arsenal, Maryland 21010	1	Director Walter Reed Army Institute of Research Walter Reed Army Medical Center Washington, D.C. 20012
1	Commanding General U.S. Army Combat Developments Command ATTN: CDCRE Fort Belvoir, Virginia 22060	1	Commanding Officer U.S. Army Medical Research Laboratory Fort Knox, Kentucky 40120
1	Commanding General U.S. Army Combat Developments Command Institute of Advanced Studies Carlisle Barracks, Pennsylvania 17013	1	Commanding Officer U.S. Naval Radiological Defense Laboratory San Francisco, California 94135
1	Commanding Officer U.S. Army Combat Developments Command Nuclear Group Fort Bliss, Texas 79916	1	Chief of Naval Research Department of the Navy Washington, D.C. 20360
1	Commandant U.S. Army Command and General Staff College Fort Leavenworth, Kansas 66027	2	Commander U.S. Naval Ordnance Laboratory White Oak Silver Spring, Maryland 20910
1	Commandant U.S. Army Air Defense School Fort Bliss, Texas 79916	1	Commander U.S. Naval Ordnance Test Station ATTN: Tech Lib China Lake, California 93557
2	Office of Civil Defense Department of the Army Washington, D.C. 20310	1	Superintendent U.S. Naval Postgraduate School ATTN: Tech Rpt Sec Monterey, California 93940
		1	Director U.S. Naval Research Laboratory Washington, D.C. 20390

# DISTRIBUTION LIST

<u>No. of</u> <u>Copies</u>	<u>Organization</u>	<u>No. of</u> <u>Copies</u>	<u>Organization</u>
1	Commandant U.S. Marine Corps School Quantico, Virginia 22134	1	Director Argonne National Laboratory 9700 South Cass Avenue Argonne, Illinois 60440
1	AFCRL (Dr. Marcus O'Day) L. G. Hanscom Fld Bedford, Massachusetts 01731	1	Director Brookhaven National Laboratory 25 Brookhaven Avenue Upton, New York 11973
1	AFWL (WLL) Kirtland AFB New Mexico 87117	1	University of California Lawrence Radiation Laboratory ATTN: Tech Lib P.O. Box 808 Livermore, California 94551
1	AUL (3T-AUL-60-118) Maxwell AFB Alabama 36112		
1	AFFDL (FDFN) Wright-Patterson AFB Ohio 45433	1	Director Los Alamos Scientific Laboratory ATTN: Tech Lib P.O. Box 1663 Los Alamos, New Mexico 87544
1	Director, Project RAND Department of the Air Force 1700 Main Street Santa Monica, California 90406	2	Director Oak Ridge National Laboratory ATTN: Dr. J. E. Strain (1 cy) Oak Ridge, Tennessee 37831
1	Director U.S. National Bureau of Standards Washington, D.C. 20234	1	Oak Ridge Institute of Nuclear Studies ATTN: Dr. H. E. Banta P.O. Box 117 Oak Ridge, Tennessee 37837
1	Director U.S. Public Health Service Washington, D.C. 20201	1	Battelle Memorial Institute 505 King Street Columbus, Ohio 43201
3	Headquarters U.S. Atomic Energy Commission ATTN: TI BM, Mr. Deal R Washington, D.C. 20545	1	General Electric Company Knolls Atomic Power Laboratory P.O. Box 1072 Schenectady, New York 13205
2	Director U.S. Atomic Energy Commission ATTN: Tech Info Div P.O. Box 62 Oak Ridge, Tennessee 37831	1	Hanford Atomic Products Operation Richland, Washington 99352



## DISTRIBUTION LIST

<u>No. of Copies</u>	<u>Organization</u>
1	Lockheed Aircraft Corporation Lockheed Nuclear Products Georgia Division Marietta, Georgia 30061
1	TRW Systems Group ATTN: Tech Lib One Space Park Redondo Beach, California 90278
1	United Nuclear Corporation White Plains, New York 10600
1	Westinghouse Electric Corporation Bettis Atomic Power Laboratory P.O. Box 355 Pittsburgh, Pennsylvania 15230
1	Massachusetts Institute of Technology 77 Massachusetts Avenue Cambridge, Massachusetts 02139

### Aberdeen Proving Ground

Ch, Tech Lib  
Air Force Ln Ofc  
Marine Corps Ln Ofc  
Navy Ln Ofc  
CDC Ln Ofc

Unclassified

## Security Classification

## DOCUMENT CONTROL DATA - R&amp;D

(Security classification of title, body of abstract and indexing annotation must be entered when the overall report is classified)

1. ORIGINATING ACTIVITY (Corporate author) U.S. Army Ballistic Research Laboratories Aberdeen Proving Ground, Maryland		2a. REPORT SECURITY CLASSIFICATION Unclassified	
		2b. GROUP	
3. REPORT TITLE  APPLICATION OF FAST NEUTRON ACTIVATION ANALYSIS TO POLYMER DOSIMETRY			
4. DESCRIPTIVE NOTES (Type of report and inclusive dates)			
5. AUTHOR(S) (Last name, first name, initial)  Coates, Arthur D.; Swisher, Joe A.			
6. REPORT DATE March 1966		7a. TOTAL NO. OF PAGES 24	7b. NO. OF REFS 11
8a. CONTRACT OR GRANT NO.		9a. ORIGINATOR'S REPORT NUMBER(S)  Report No. 1312	
b. PROJECT NO. RDT&E 1N022601A086			
c.		9b. OTHER REPORT NO(S) (Any other numbers that may be assigned this report)	
d.			
10. AVAILABILITY/LIMITATION NOTICES  Distribution of this document is unlimited.			
11. SUPPLEMENTARY NOTES		12. SPONSORING MILITARY ACTIVITY U.S. Army Materiel Command Washington, D.C.	
13. ABSTRACT  Fast neutron activation analysis was used to measure the gamma dose absorbed by a polymer. A series of low density polyethylene films was irradiated by cobalt-60 gamma rays in the presence of air. The gamma doses ranged from $1 \times 10^5$ to $4.7 \times 10^7$ rads. The oxygen content of the polymer films increased with increasing dose. The increase in oxygen content of each film was monitored at $1720 \text{ cm}^{-1}$ by infrared spectroscopy. The oxygen increase also was monitored by measuring the activation of oxygen-16. This latter method employs fast neutron activation of oxygen-16 and the measurement of the emitted radiation to determine gamma dose. The results of both methods of dosimetry are presented, and the improvements offered by the neutron activation technique are indicated.			

14. KEY WORDS	LINK A		LINK B		LINK C	
	ROLE	WT	ROLE	WT	ROLE	WT
Gamma Dosimetry Polymer Dosimetry Fast Neutron Activation Oxygen Analysis						

#### INSTRUCTIONS

1. **ORIGINATING ACTIVITY:** Enter the name and address of the contractor, subcontractor, grantee, Department of Defense activity or other organization (*corporate author*) issuing the report.

2a. **REPORT SECURITY CLASSIFICATION:** Enter the overall security classification of the report. Indicate whether "Restricted Data" is included. Marking is to be in accordance with appropriate security regulations.

2b. **GROUP:** Automatic downgrading is specified in DoD Directive 5200.10 and Armed Forces Industrial Manual. Enter the group number. Also, when applicable, show that optional markings have been used for Group 3 and Group 4 as authorized.

3. **REPORT TITLE:** Enter the complete report title in all capital letters. Titles in all cases should be unclassified. If a meaningful title cannot be selected without classification, show title classification in all capitals in parenthesis immediately following the title.

4. **DESCRIPTIVE NOTES:** If appropriate, enter the type of report, e.g., interim, progress, summary, annual, or final. Give the inclusive dates when a specific reporting period is covered.

5. **AUTHOR(S):** Enter the name(s) of author(s) as shown on or in the report. Enter last name, first name, middle initial. If military, show rank and branch of service. The name of the principal author is an absolute minimum requirement.

6. **REPORT DATE:** Enter the date of the report as day, month, year; or month, year. If more than one date appears on the report, use date of publication.

7a. **TOTAL NUMBER OF PAGES:** The total page count should follow normal pagination procedures, i.e., enter the number of pages containing information.

7b. **NUMBER OF REFERENCES:** Enter the total number of references cited in the report.

8a. **CONTRACT OR GRANT NUMBER:** If appropriate, enter the applicable number of the contract or grant under which the report was written.

8b, 8c, & 8d. **PROJECT NUMBER:** Enter the appropriate military department identification, such as project number, subproject number, system numbers, task number, etc.

9a. **ORIGINATOR'S REPORT NUMBER(S):** Enter the official report number by which the document will be identified and controlled by the originating activity. This number must be unique to this report.

9b. **OTHER REPORT NUMBER(S):** If the report has been assigned any other report numbers (*either by the originator or by the sponsor*), also enter this number(s).

10. **AVAILABILITY/LIMITATION NOTICES:** Enter any limitations on further dissemination of the report, other than those imposed by security classification, using standard statements such as:

- (1) "Qualified requesters may obtain copies of this report from DDC."
- (2) "Foreign announcement and dissemination of this report by DDC is not authorized."
- (3) "U. S. Government agencies may obtain copies of this report directly from DDC. Other qualified DDC users shall request through \_\_\_\_\_."
- (4) "U. S. military agencies may obtain copies of this report directly from DDC. Other qualified users shall request through \_\_\_\_\_."
- (5) "All distribution of this report is controlled. Qualified DDC users shall request through \_\_\_\_\_."

If the report has been furnished to the Office of Technical Services, Department of Commerce, for sale to the public, indicate this fact and enter the price, if known.

11. **SUPPLEMENTARY NOTES:** Use for additional explanatory notes.

12. **SPONSORING MILITARY ACTIVITY:** Enter the name of the departmental project office or laboratory sponsoring (*paying for*) the research and development. Include address.

13. **ABSTRACT:** Enter an abstract giving a brief and factual summary of the document indicative of the report, even though it may also appear elsewhere in the body of the technical report. If additional space is required, a continuation sheet shall be attached.

It is highly desirable that the abstract of classified reports be unclassified. Each paragraph of the abstract shall end with an indication of the military security classification of the information in the paragraph, represented as (TS), (S), (C), or (U).

There is no limitation on the length of the abstract. However, the suggested length is from 150 to 225 words.

14. **KEY WORDS:** Key words are technically meaningful terms or short phrases that characterize a report and may be used as index entries for cataloging the report. Key words must be selected so that no security classification is required. Identifiers, such as equipment model designation, trade name, military project code name, geographic location, may be used as key words but will be followed by an indication of technical context. The assignment of links, rules, and weights is optional.